
# NeoMed NEOCONNECT™

## NeoConnect™ Oral/Enteral Syringe with ENFit® connector For Pharmacy Use; Instructions for Use

- CZ:** Orální / enterální stříkačka NeoConnect s konektorem ENFit; Pro lékárenské použití ; návod k použití  
**DK:** NeoConnect oral/enteral sprøjte med ENFit-stik; Til apoteksbrug; Brugershåndbog  
**NL:** NeoConnect orale/enterale spuit met ENFit-connector; Voor apotheekgebruik; gebruiksaanwijzing  
**FI:** NeoConnect- oraallinen/enteraalinen ruisku ENFit-liittimellä; Farmasiakäyttöön; Käyttöohjeet  
**FR:** Seringue orale/entérale NeoConnect avec connecteur ENFit ; Pour un usage pharmaceutique ; mode d'emploi  
**DE:** NeoConnect orale/enterale Spritze mit ENFit-Konnektor; Für die pharmazeutische Verwendung; Gebrauchsanweisung  
**IT:** Siringa orale/enterale NeoConnect con connettore ENFit; Per uso farmaceutico; Istruzioni per l'uso  
**NO:** NeoConnect oral/enteral sprøyte med ENFit-kobling. For apotekbruk. Bruksanvisning.  
**PL:** Strzykawka doustna/dojelitowa NeoConnect ze złączem ENFit; Do użytku aptecznego; Instrukcja użytkowania  
**PT:** Seringa oral/etérica NeoConnect com conector ENFit; Para utilização farmacêutica; Instruções de utilização  
**RU:** Оральный/энтеральный шприц NeoConnect с разъемом ENFit, Для применения в фармацевтике; инструкция по использованию  
**ES:** Jeringa oral/enteral NeoConnect con conector ENFit; Para uso en farmacia; Instrucciones de uso  
**SE:** NeoConnect oral/enteral spruta med ENFit-anslutning, För användning i apotek, bruksanvisning.

Figure 1; Obrázek 1; Figur 1; Figura 1; Kuva 1; Obrázek 1  
Afbeelding 1; Rysunek 1; Рис. 1; Figure 1

### Indications for Use:

The device is indicated for use as a dispenser, a measuring device and a fluid transfer device. It is used to deliver fluids into the body orally or enterally. It is intended to be used in clinical or home care settings by users ranging from clinicians to laypersons (under the supervision of a clinician) in all age groups.

### Contraindications:

- The NeoMed NeoConnect Syringe is intended to connect to products with ENFit compatible male connectors.
- This device is intended for oral / enteral applications only. Do not use this device with small bore connectors from other healthcare applications as it has the potential to misconnect.

### Warnings:

- For syringe sizes 6mL and below, to ensure dose accuracy, the moat (See Figure 1) must be free of fluids prior to patient administration.
- The syringe is intended for one-time use. Each syringe should be used for one dispensing on one patient.
- Do not autoclave or re-use syringe as this may increase the risk of leakage or biological infection.
- Tip caps are a choking hazard. Keep out of reach of children.
- Do not modify device as this may lead to leakage, inadequate nutrition delivery or patient harm.
- Do not use device if syringe tip is visibly blocked as this may lead to inadequate nutrition delivery.
- Do not use if device is damaged as this may lead to device contamination or inadequate nutrition delivery.

### Directions:

- Follow physician's instructions or hospital protocol in opening package and use of contents.
- Always use technique that minimizes risk of infection. Syringe should be connected to a male ENFit compatible connector during filling and removal of air bubbles to reduce the risk of dosing volume inaccuracies.
- Fully depress the plunger before filling syringe. Draw liquid (e.g. formula or mother's milk) into syringe. For accurate measurement, use the black leading edge of the gasket toward the plunger tip. Do not fill and attempt delivery past the maximum graduation.
- To expel excess air bubbles, hold syringe with tip pointing up. Tap side of syringe to allow bubbles to rise toward tip. Draw back slightly on plunger, then push plunger upward to eject air in barrel and tip. Repeat until bubbles are eliminated and re-verify dose.
- For syringe sizes 6ml and below, ensure the moat (See Figure 1) is free of fluids and re-verify dose (repeat as necessary).
- Secure the tip cap on the syringe and label until ready to use (label not included).
- Before feeding, to minimize potential for clogging, mix any supplements or additives thoroughly to achieve uniform consistency.
- For direct oral administration, DoseMate™ may be attached to syringe (sold separately) for patient comfort.
- Dispose of syringe as required by hospital protocols and local or national regulations.

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DK

Manufactured by NeoMed Incorporated,  
100 Londonderry Ct, Suite 112, Woodstock, GA 30188 USA  
(888) 876-2225; www.neomedinc.com

Emergo Europe, Molenstraat 15, 2513 BH,  
The Hague, The Netherlands

Australian Sponsor: Emergo Australia  
Level 20, Tower II Darling Park  
201 Sussex Street, Sydney, NSW 2000  
Australia

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