

Medi-Dose / EPS


Unit-Dose, Bar Coding, Pharmacy and Nursing Supply Experts!

Neuromuscular Blocking Agents Warning Label Image Install


These images © 2018 Medi-Dose®, Inc. / EPS®, Inc.
Expressed permissions are granted to use and print these images with MILT by Medi-Dose Software only.
Any other reproduction or distribution of these images is forbidden.

Step 1: [DOWNLOAD HERE](#), open the file and press OK to agree with the terms

Step 2: CHOOSE UNZIP



Step 3: PRESS OK




Step 4: Open MILT 4


Step 5: Select “3. Template” and at the bottom of the screen and select the IMAGE button




Step 6: Choose IMPORT



Step 7: Browse to C:\NBAW\



Step 8: Select the image to import then choose OPEN



The image is now part of your library.

Simply select that image to insert it into your template, resize and relocate it, then save the template.

Medi-Dose® / EPS®
70 Industrial Drive Ivyland, PA 18974
MediDose.com 800.523.8966