

FOR IMMEDIATE RELEASE:

CONTACT:

Robert Braverman
Medi-Dose, Inc. / EPS, Inc.
800-523-8966 – Phone
800-323-8966 – Fax
rbraverman@medidose.com
www.medidose.com

Medi-Dose® / EPS® Announces New ENFit® at Home Syringes

Ivyland, PA, Spring 2020 – The cost and accessibility of ENFit syringes for discharged patients has challenged the overall adoption of this new standard for enteral feeding...until now. To address these concerns, Medi-Dose/EPS is pleased to announce their new **ENFit at Home Syringes**, extending its existing line of ENFit Compliant and Compatible Products.

These new **ENFit at Home Syringes** are approved to be hand washed and reused up to 20 times, providing a cost-effective and eco-friendly alternative to single-use ENFit syringes. Smaller case quantities make them easier for individual families to afford. Since a prescription is not required, they are also easier to obtain.

All 4 sizes (3, 6, 12 and 60 mL) feature a durable polypropylene plunger with a silicone O-Ring gasket, a clear barrel and black, easy-to read, metric only graduations. These graduations comply with ISMP Best Standards Guidelines.

“Reducing misconnections with medical device tubing is vital to improving patient safety. Our new ENFit at Home syringes help families deliver enteral medication safely, securely and economically,” says Bob Braverman, President.

For more information on the Medi-Dose/EPS line of ENFit at Home Syringes and other NeoConnect ENFit products, contact Robert Braverman or visit www.medidose.com.

ABOUT MEDI-DOSE, INC. / EPS, INC. – For 50 years, Medi-Dose / EPS has been working with pharmacists and health care professionals to develop cost-effective solutions that promote medication safety and error prevention.

- END -