

FOR IMMEDIATE RELEASE:

Medi-Dose® , Inc. / EPS® , Inc.
Now Markets Steri-Tamp® Tamper-Clear™ Syringe Seals

Ivyland, PA, Fall 2017 – To complement their growing lines of IV and tamper-evident products, Medi-Dose/EPS is pleased to add the new Steri-Tamp Tamper-Clear Syringe Seal.

The Tamper-Clear Syringe Seal provides a tamper-evident seal with clear visibility of syringe graduations. It's easily applied over the syringe cap and maintains the integrity of the medication and packaging...even when refrigerated. Once removed, the seal can't be reapplied. Removing the syringe cap destroys the seal.

Designed to withstand freezing temperatures down to -20 degrees Centigrade, Steri-Tamp®'s powerful adhesive seal has been shown to be over three times stronger than others on the market. It provides a 100% sterile barrier. All Steri-Tamp® products are packaged in a convenient USP <797> compliant dispenser.

"We're always looking for ways to help pharmacists and nurses streamline their practices but maintain control and protection over the meds they dispense," says Bob Braverman, Director of Marketing. "The new Tamper-Clear Syringe Seal as well as the rest of Steri-Tamp® line of bag, vial and syringe seals help practitioners accomplish this in a cost-effective manner. These products are an excellent addition to the fight against contamination, tampering and accidental double dosing."

For more information on the EPS Steri-Tamp® line of tamper-evident IV seals, please visit www.MediDose.com, call 800-523-8966 or e-mail info@medidose.com.

ABOUT MEDI-DOSE, INC. / EPS, INC. – For over 45 years, Medi-Dose / EPS has been working with pharmacists to develop cost-effective solutions that promote medication safety and error prevention.

###