


Medi-Dose[®] EPS[®]

Unit Dose,
Bar Coding,
Pharmacy &
Nursing Supply
Experts!

STERI-TAMP[®] IV Seals



The Only Tamper-Evident Seal that Provides a 100% Sterile Barrier

- Protects IV Bags, Vials and Syringes from Contamination and Double Dosing
- Removing the Seal Leaves Printed "OPENED" Warning
- Foil Outer Layer Shows Needle Punctures
- Maintains Adhesion When Refrigerated or Frozen (Down to -20 Degrees C)
- IV Bags with Unbroken Seals May be Returned and Reused
- Packaged in Easy to Use Dispensers with 1,000 Seals per Roll

TAMPER-CLEAR SYRINGE SEAL[™]



Same Great Tamper-Evidence with a See-Through Seal!

MediDose.com

800.523.8966

Exceptional Products & Service

Unit Dose,
Bar Coding,
Pharmacy &
Nursing Supply
Experts!

Medi-Dose System



Unit Dose



Liquid Packaging



Steri-Dropper Bottles



TampAlert[™] System



Bags



IV Pharmacy/Nursing



Labels



Tamper Evidence



Storage & Transport



ENFit[™] Products



Cleanroom Supplies



Patient Care



Laboratory



Identification Products



Software

